


INSTRUCTIONS

If you have a letter inviting you to do the test with 2 labels, follow steps 1 and 2, then go to step 3.
If you do not have this letter, go straight to step 1bis below.




1 Stick the large label on the identification form (located in the central part of the kit) and **complete the form**.



2 Record the test date on the small label and **stick it** on the flat part of the tube where "name" and "date" appear.



1 bis Carefully complete the identification form, including the label in the central part of the kit.




2 bis Stick the label, on the flat part of the tube where "name" and "date" appear. Proceed to step 3.



3 Stick the adhesive paper for collecting the stools on the toilet seat.
Pat it in the middle to create a dip.


IMPORTANT: for the test to work, there must be no contact with any liquid, so no urine, no bleach, nothing.



4 Open the tube by unscrewing the top.



5 Use the stick to scratch the surface of the stool in several places.



6 The end of the stick has ridges. The ridges must be completely covered in stool.



7 Close the tube tightly and give it a shake.



8 Check that you have stuck the completed label on the tube.
Put the tube into its protective sleeve.



9 Take the completed identification form and put everything in the prepaid envelope.
Close the envelope.



10 Post the envelope the following day at the latest.

You should receive your results by post within two weeks.

For online results, please register at the following website: www.resultat-depistage.fr

FOR FURTHER INFORMATION?

- Ask your GP, chemist or nurse.
- On line at: jefaismondepistage.fr
- Call:
0 805 123 124 (local call cost)
Monday-Friday, 9:00 am-7:00 pm
and Saturday, 9:00 am-2:00 pm.
- Contact the department responsible for screening in your area.

FAQ

- **Within what period of time must the test be used?**
The expiry date is on the tube.
- **What if the tube gets lost or damaged?**
Ask your GP, chemist or nurse for a new test.
- **Something went wrong with the sample collection, what should I do?**
Contact your GP, chemist or nurse; he will advise you accordingly and give you a new test.